

Pig CSTB/Cystatin-B ELISA Kit

Catalog No: #EK19824



Orders: order@signalwayantibody.com
Support: tech@signalwayantibody.com

Description

Product Name	Pig CSTB/Cystatin-B ELISA Kit
Brief Description	ELISA Kits
Applications	ELISA
Species Reactivity	Pig
Specificity	Natural and recombinant Pig Cystatin-B
Target Name	Cystatin-B
Other Names	Stefin-B,CST6,STFB
Accession No.	Q29290
Uniprot	Q29290

Application Details

Assay Type: SandwichDetect Range:0.312-20ng/mL Sensitivity:0.156ng/mL Sample Type:Serum, plasma, tissue homogenates and other biological fluids.
Sample Volume:100ul Assay Time:3-5h Detection wavelength:450 nmgene_full_name:CSTB

Note: This product is for in vitro research use only