

Rat Cstb/Cystatin-B ELISA Kit

Catalog No: #EK19825

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Rat Cstb/Cystatin-B ELISA Kit
Brief Description	ELISA Kits
Applications	ELISA
Species Reactivity	Rt
Specificity	Natural and recombinant Rat Cystatin-B
Target Name	Cystatin-B
Other Names	Cystatin-beta,Liver thiol proteinase inhibitor,Stefin-B,Cst6
Accession No.	P01041
Uniprot	P01041
GeneID	25308;

Application Details

Assay Type: Sandwich Detect Range:0.312-20ng/mL Sensitivity:0.156ng/mL Sample Type:Serum, plasma, tissue homogenates and other biological fluids.

Sample Volume:100ul Assay Time:3-5h Detection wavelength:450 nm gene_full_name:

Note: This product is for in vitro research use only