

Human CALU/Calumenin ELISA Kit

Catalog No: #EK19826

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Human CALU/Calumenin ELISA Kit
Brief Description	ELISA Kits
Applications	ELISA
Species Reactivity	Hu
Specificity	Natural and recombinant Human Calumenin
Target Name	Calumenin
Other Names	Crocalbin,IEF SSP 9302
Accession No.	O43852
Uniprot	O43852
GeneID	813;

Application Details

Assay Type: Sandwich Detect Range:0.312-20ng/mL Sensitivity:0.156ng/mL Sample Type:Serum, plasma, tissue homogenates and other biological fluids.

Sample Volume:100ul Assay Time:3-5h Detection wavelength:450 nm gene_full_name:CALU

Note: This product is for in vitro research use only