

Mouse Calu/Calumenin ELISA Kit

Catalog No: #EK19827

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Mouse Calu/Calumenin ELISA Kit
Brief Description	ELISA Kits
Applications	ELISA
Species Reactivity	Ms
Specificity	Natural and recombinant Mouse Calumenin
Target Name	Calumenin
Other Names	Crocalbin
Accession No.	O35887
Uniprot	O35887
GeneID	12321;

Application Details

Assay Type: Sandwich Detect Range:0.312-20ng/mL Sensitivity:0.156ng/mL Sample Type:Serum, plasma, tissue homogenates and other biological fluids.

Sample Volume:100ul Assay Time:3-5h Detection wavelength:450 nm gene_full_name:Calu

Note: This product is for in vitro research use only