

Rat Calu/Calumenin ELISA Kit

Catalog No: #EK19828

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Rat Calu/Calumenin ELISA Kit
Brief Description	ELISA Kits
Applications	ELISA
Species Reactivity	Rt
Specificity	Natural and recombinant Rat Calumenin
Target Name	Calumenin
Other Names	CBP-50, Crocalbin
Accession No.	O35783
Uniprot	O35783

Application Details

Assay Type: Sandwich Detect Range: 1.56-100ng/mL Sensitivity: 0.88ng/mL Sample Type: Serum, plasma, tissue homogenates and other biological fluids.
Sample Volume: 100ul Assay Time: 3-5h Detection wavelength: 450 nm gene_full_name: Calu

Note: This product is for in vitro research use only