

Mouse Cfl2/Cofilin-2 ELISA Kit

Catalog No: #EK19830

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Mouse Cfl2/Cofilin-2 ELISA Kit
Brief Description	ELISA Kits
Applications	ELISA
Species Reactivity	Ms
Specificity	Natural and recombinant Mouse Cofilin-2
Target Name	Cofilin-2
Other Names	Cofilin, muscle isoform
Accession No.	P45591
Uniprot	P45591
GeneID	12632;

Application Details

Assay Type: Sandwich Detect Range: 1.25-80ng/mL Sensitivity: 0.633ng/mL Sample Type: Serum, plasma, tissue homogenates and other biological fluids.
Sample Volume: 100ul Assay Time: 3-5h Detection wavelength: 450 nm gene_full_name: Cfl2

Note: This product is for in vitro research use only