

Mouse Cd3e/T-cell surface glycoprotein CD3 epsilon chain ELISA Kit

Catalog No: #EK19831

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Mouse Cd3e/T-cell surface glycoprotein CD3 epsilon chain ELISA Kit
Brief Description	ELISA Kits
Applications	ELISA
Species Reactivity	Ms
Specificity	Natural and recombinant Mouse T-cell surface glycoprotein CD3 epsilon chain
Target Name	T-cell surface glycoprotein CD3 epsilon chain
Other Names	T-cell surface antigen T3/Leu-4 epsilon chain,CD3e
Accession No.	P22646
Uniprot	P22646
GeneID	12501;

Application Details

Assay Type: Sandwich Detect Range:0.312-20ng/mL Sensitivity:0.135ng/mL Sample Type:Serum, plasma, tissue homogenates and other biological fluids.

Sample Volume:100ul Assay Time:3-5h Detection wavelength:450 nm gene_full_name: Cd3e

Note: This product is for in vitro research use only