

Human CD19/B-lymphocyte antigen CD19 ELISA Kit

Catalog No: #EK19832



Orders: order@signalwayantibody.com
Support: tech@signalwayantibody.com

Description

Product Name	Human CD19/B-lymphocyte antigen CD19 ELISA Kit
Brief Description	ELISA Kits
Applications	ELISA
Species Reactivity	Hu
Specificity	Natural and recombinant Human B-lymphocyte antigen CD19
Target Name	B-lymphocyte antigen CD19
Other Names	B-lymphocyte surface antigen B4,Differentiation antigen CD19,T-cell surface antigen Leu-12,CD19
Accession No.	P15391
Uniprot	P15391
GeneID	930;

Application Details

Assay Type: SandwichDetect Range:0.312-20ng/mL Sensitivity:0.168ng/mL Sample Type:Serum, plasma, tissue homogenates and other biological fluids.
Sample Volume:100ul Assay Time:3-5h Detection wavelength:450 nmgene_full_name:CD19

Note: This product is for in vitro research use only