

Mouse Cd19/B-lymphocyte antigen CD19 ELISA Kit

Catalog No: #EK19833

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Mouse Cd19/B-lymphocyte antigen CD19 ELISA Kit
Brief Description	ELISA Kits
Applications	ELISA
Species Reactivity	Ms
Specificity	Natural and recombinant Mouse B-lymphocyte antigen CD19
Target Name	B-lymphocyte antigen CD19
Other Names	Differentiation antigen CD19,CD19
Accession No.	P25918
Uniprot	P25918
GeneID	12478;

Application Details

Assay Type: Sandwich Detect Range:0.312-20ng/ml Sensitivity:0.171ng/mL Sample Type:Serum, plasma, tissue homogenates and other biological fluids.
Sample Volume:100ul Assay Time:3-5h Detection wavelength:450 nm gene_full_name: Cd19

Note: This product is for in vitro research use only