

Human CDK7/Cyclin-dependent kinase 7 ELISA Kit

Catalog No: #EK19834

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Human CDK7/Cyclin-dependent kinase 7 ELISA Kit
Brief Description	ELISA Kits
Applications	ELISA
Species Reactivity	Hu
Specificity	Natural and recombinant Human Cyclin-dependent kinase 7
Target Name	Cyclin-dependent kinase 7
Other Names	39 kDa protein kinase,p39 Mo15,CDK-activating kinase 1,Cell division protein kinase 7,Serine/threonine-protein kinase 1,TFIIH basal transcription factor complex kinase subunit,CAK,CAK1,CDKN7,MO15,STK1
Accession No.	P50613
Uniprot	P50613
GeneID	1022;

Application Details

Assay Type: Sandwich Detect Range:0.312-20ng/mL Sensitivity:0.174ng/mL Sample Type:Serum, plasma, tissue homogenates and other biological fluids.
Sample Volume:100ul Assay Time:3-5h Detection wavelength:450 nm gene_full_name:CDK7

Note: This product is for in vitro research use only