

## Rat Cdk7/Cyclin-dependent kinase 7 ELISA Kit

Catalog No: #EK19836

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

## Description

Product Name	Rat Cdk7/Cyclin-dependent kinase 7 ELISA Kit
Brief Description	ELISA Kits
Applications	ELISA
Species Reactivity	Rt
Specificity	Natural and recombinant Rat Cyclin-dependent kinase 7
Target Name	Cyclin-dependent kinase 7
Other Names	39 protein kinase,P39 Mo15,CDK-activating kinase 1,Cell division protein kinase 7,TFIIH basal transcription factor complex kinase subunit,Cak,Cak1,Mo15
Accession No.	P51952
Uniprot	P51952

## Application Details

Assay Type: Sandwich Detect Range:78-5000pg/mL Sensitivity:28pg/mL Sample Type:Serum, plasma, tissue homogenates and other biological fluids. Sample Volume:100ul Assay Time:3-5h Detection wavelength:450 nm gene\_full\_name:Cdk7

Note: This product is for in vitro research use only