

## Chicken CDK1/Cyclin-dependent kinase 1 ELISA Kit

Catalog No: #EK19839

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

## Description

Product Name	Chicken CDK1/Cyclin-dependent kinase 1 ELISA Kit
Brief Description	ELISA Kits
Applications	ELISA
Species Reactivity	Chicken
Specificity	Natural and recombinant Chicken Cyclin-dependent kinase 1
Target Name	Cyclin-dependent kinase 1
Other Names	CDK1, Cell division control protein 2 homolog, Cell division protein kinase 1, p34 protein kinase, CDC2
Accession No.	P13863
Uniprot	P13863
GeneID	396252;

## Application Details

Assay Type: Sandwich Detect Range: 0.312-20ng/ml Sensitivity: 0.18ng/mL Sample Type: Serum, plasma, tissue homogenates and other biological fluids.  
Sample Volume: 100ul Assay Time: 3-5h Detection wavelength: 450 nm gene\_full\_name: CDK1

Note: This product is for in vitro research use only