

## Rat Cdk1/Cyclin-dependent kinase 1 ELISA Kit

Catalog No: #EK19842

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

## Description

Product Name	Rat Cdk1/Cyclin-dependent kinase 1 ELISA Kit
Brief Description	ELISA Kits
Applications	ELISA
Species Reactivity	Rt
Specificity	Natural and recombinant Rat Cyclin-dependent kinase 1
Target Name	Cyclin-dependent kinase 1
Other Names	CDK1, Cell division control protein 2 homolog, Cell division protein kinase 1, p34 protein kinase, Cdc2, Cdc2a, Cdkn1
Accession No.	P39951
Uniprot	P39951
GeneID	54237;

## Application Details

Assay Type: Sandwich Detect Range: 0.78-50ng/mL Sensitivity: 0.396ng/mL Sample Type: Serum, plasma, tissue homogenates and other biological fluids.  
Sample Volume: 100ul Assay Time: 3-5h Detection wavelength: 450 nm gene\_full\_name: Cdk1

Note: This product is for in vitro research use only