

## Bovine CDK2/Cyclin-dependent kinase 2 ELISA Kit

Catalog No: #EK19843

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

## Description

Product Name	Bovine CDK2/Cyclin-dependent kinase 2 ELISA Kit
Brief Description	ELISA Kits
Applications	ELISA
Species Reactivity	Bovine
Specificity	Natural and recombinant Bovine Cyclin-dependent kinase 2
Target Name	Cyclin-dependent kinase 2
Other Names	Cell division protein kinase 2,CDKN2
Accession No.	Q5E9Y0
Uniprot	Q5E9Y0
GeneID	519217;

## Application Details

Assay Type: Sandwich Detect Range: 78-5000pg/mL Sensitivity: 39pg/mL Sample Type: Serum, plasma, tissue homogenates and other biological fluids. Sample Volume: 100ul Assay Time: 3-5h Detection wavelength: 450 nm gene\_full\_name: CDK2

Note: This product is for in vitro research use only