

Rat Cdk2/Cyclin-dependent kinase 2 ELISA Kit

Catalog No: #EK19844

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Rat Cdk2/Cyclin-dependent kinase 2 ELISA Kit
Brief Description	ELISA Kits
Applications	ELISA
Species Reactivity	Rt
Specificity	Natural and recombinant Rat Cyclin-dependent kinase 2
Target Name	Cyclin-dependent kinase 2
Other Names	Cell division protein kinase 2,Cdkn2
Accession No.	Q63699
Uniprot	Q63699

Application Details

Assay Type: Sandwich Detect Range: 78-5000pg/mL Sensitivity: 39pg/mL Sample Type: Serum, plasma, tissue homogenates and other biological fluids. Sample Volume: 100ul Assay Time: 3-5h Detection wavelength: 450 nm gene_full_name: Cdk2

Note: This product is for in vitro research use only