

Human ACVRL1/Serine/threonine-protein kinase receptor R3 ELISA Kit

Catalog No: #EK19846

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Human ACVRL1/Serine/threonine-protein kinase receptor R3 ELISA Kit
Brief Description	ELISA Kits
Applications	ELISA
Species Reactivity	Hu
Specificity	Natural and recombinant Human Serine/threonine-protein kinase receptor R3
Target Name	Serine/threonine-protein kinase receptor R3
Other Names	SKR3,Activin receptor-like kinase 1,ALK-1,TGF-B superfamily receptor type I,TSR-I,ACVRLK1,ALK1
Accession No.	P37023
Uniprot	P37023
GeneID	94;

Application Details

Assay Type: Sandwich Detect Range:0.78-50ng/mL Sensitivity:0.397ng/mL Sample Type:Serum, plasma, tissue homogenates and other biological fluids.
Sample Volume:100ul Assay Time:3-5h Detection wavelength:450 nm gene_full_name:ACVRL1

Note: This product is for in vitro research use only