

Mouse Acvr11/Serine/threonine-protein kinase receptor R3 ELISA Kit



Catalog No: #EK19847

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Mouse Acvr11/Serine/threonine-protein kinase receptor R3 ELISA Kit
Brief Description	ELISA Kits
Applications	ELISA
Species Reactivity	Ms
Specificity	Natural and recombinant Mouse Serine/threonine-protein kinase receptor R3
Target Name	Serine/threonine-protein kinase receptor R3
Other Names	SKR3,Activin receptor-like kinase 1,ALK-1,TGF-B superfamily receptor type I,TSR-I,Acvr1k1,Alk-1
Accession No.	Q61288
Uniprot	Q61288
GeneID	11482;

Application Details

Assay Type: SandwichDetect Range:0.312-20ng/mL Sensitivity:0.098ng/mL Sample Type:Serum, plasma, tissue homogenates and other biological fluids.

Sample Volume:100ul Assay Time:3-5h Detection wavelength:450 nmgene_full_name:Acvr11

Note: This product is for in vitro research use only