

Rat Acvr11/Serine/threonine-protein kinase receptor R3 ELISA Kit



Catalog No: #EK19848

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

| | |
|--------------------|---|
| Product Name | Rat Acvr11/Serine/threonine-protein kinase receptor R3 ELISA Kit |
| Brief Description | ELISA Kits |
| Applications | ELISA |
| Species Reactivity | Rt |
| Specificity | Natural and recombinant Rat Serine/threonine-protein kinase receptor R3 |
| Target Name | Serine/threonine-protein kinase receptor R3 |
| Other Names | SKR3,TGF-B superfamily receptor type I,TSR-I,Acvr1k1 |
| Accession No. | P80203 |
| Uniprot | P80203 |
| GeneID | 25237; |

Application Details

Assay Type: Sandwich Detect Range:0.78-50ng/mL Sensitivity:0.398ng/mL Sample Type:Serum, plasma, tissue homogenates and other biological fluids.
Sample Volume:100ul Assay Time:3-5h Detection wavelength:450 nm gene_full_name:Acvr11

Note: This product is for in vitro research use only