

## Rat Acan/Aggrecan core protein ELISA Kit

Catalog No: #EK19851

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

## Description

Product Name	Rat Acan/Aggrecan core protein ELISA Kit
Brief Description	ELISA Kits
Applications	ELISA
Species Reactivity	Rt
Specificity	Natural and recombinant Rat Aggrecan core protein
Target Name	Aggrecan core protein
Other Names	Cartilage-specific proteoglycan core protein,CSPCP,Agc,Agc1
Accession No.	P07897
Uniprot	P07897
GeneID	58968;

## Application Details

Assay Type: Sandwich Detect Range:0.156-10ng/ml Sensitivity:0.091ng/mL Sample Type:Serum, plasma, tissue homogenates and other biological fluids.  
Sample Volume:100ul Assay Time:3-5h Detection wavelength:450 nm gene\_full\_name:Acan

Note: This product is for in vitro research use only