

## Mouse Bbc3/Bcl-2-binding component 3 ELISA Kit

Catalog No: #EK19853

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

## Description

Product Name	Mouse Bbc3/Bcl-2-binding component 3 ELISA Kit
Brief Description	ELISA Kits
Applications	ELISA
Species Reactivity	Ms
Specificity	Natural and recombinant Mouse Bcl-2-binding component 3
Target Name	Bcl-2-binding component 3
Other Names	p53 up-regulated modulator of apoptosis,Puma
Accession No.	Q99ML1
Uniprot	Q99ML1
GeneID	170770;

## Application Details

Assay Type: Sandwich Detect Range:0.312-20ng/mL Sensitivity:0.165ng/mL Sample Type:Serum, plasma, tissue homogenates and other biological fluids.

Sample Volume:100ul Assay Time:3-5h Detection wavelength:450 nm gene\_full\_name:Bbc3

Note: This product is for in vitro research use only