

Rat Bbc3/Bcl-2-binding component 3 ELISA Kit

Catalog No: #EK19854

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Rat Bbc3/Bcl-2-binding component 3 ELISA Kit
Brief Description	ELISA Kits
Applications	ELISA
Species Reactivity	Rt
Specificity	Natural and recombinant Rat Bcl-2-binding component 3
Target Name	Bcl-2-binding component 3
Other Names	p53 up-regulated modulator of apoptosis,Puma
Accession No.	Q80ZG6
Uniprot	Q80ZG6
GeneID	317673;

Application Details

Assay Type: Sandwich Detect Range:78-5000pg/mL Sensitivity:40pg/mL Sample Type:Serum, plasma, tissue homogenates and other biological fluids. Sample Volume:100ul Assay Time:3-5h Detection wavelength:450 nm gene_full_name:Bbc3

Note: This product is for in vitro research use only