

# Human CHRNA4/Neuronal acetylcholine receptor subunit alpha-4 ELISA Kit



Catalog No: #EK19855

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

## Description

Product Name	Human CHRNA4/Neuronal acetylcholine receptor subunit alpha-4 ELISA Kit
Brief Description	ELISA Kits
Applications	ELISA
Species Reactivity	Hu
Specificity	Natural and recombinant Human Neuronal acetylcholine receptor subunit alpha-4
Target Name	Neuronal acetylcholine receptor subunit alpha-4
Other Names	NACRA4
Accession No.	P43681
Uniprot	P43681
GeneID	1137;

## Application Details

Assay Type: Sandwich Detect Range: 0.625-40ng/mL Sensitivity: 0.341ng/mL Sample Type: Serum, plasma, tissue homogenates and other biological fluids.

Sample Volume: 100ul Assay Time: 3-5h Detection wavelength: 450 nm gene\_full\_name: CHRNA4

Note: This product is for in vitro research use only