

Mouse Angptl2/Angiopoietin-related protein 2 ELISA Kit



Catalog No: #EK19856

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Mouse Angptl2/Angiopoietin-related protein 2 ELISA Kit
Brief Description	ELISA Kits
Applications	ELISA
Species Reactivity	Ms
Specificity	Natural and recombinant Mouse Angiopoietin-related protein 2
Target Name	Angiopoietin-related protein 2
Other Names	Angiopoietin-like protein 2,Arp2
Accession No.	Q9R045
Uniprot	Q9R045
GeneID	26360;

Application Details

Assay Type: SandwichDetect Range:78-5000pg/mL Sensitivity:36pg/mL Sample Type:Serum, plasma, tissue homogenates and other biological fluids. Sample Volume:100ul Assay Time:3-5h Detection wavelength:450 nmgene_full_name:Angptl2

Note: This product is for in vitro research use only