

## Mouse Ccr4/C-C chemokine receptor type 4 ELISA Kit

Catalog No: #EK19857

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

## Description

Product Name	Mouse Ccr4/C-C chemokine receptor type 4 ELISA Kit
Brief Description	ELISA Kits
Applications	ELISA
Species Reactivity	Ms
Specificity	Natural and recombinant Mouse C-C chemokine receptor type 4
Target Name	C-C chemokine receptor type 4
Other Names	C-C CKR-4, CD194, Cmkbr4
Accession No.	P51680
Uniprot	P51680
GeneID	12773;

## Application Details

Assay Type: Sandwich Detect Range: 78-5000pg/mL Sensitivity: 34pg/mL Sample Type: Serum, plasma, tissue homogenates and other biological fluids. Sample Volume: 100ul Assay Time: 3-5h Detection wavelength: 450 nm gene\_full\_name: Ccr4

Note: This product is for in vitro research use only