

# Human CHRNA2/Neuronal acetylcholine receptor subunit alpha-2 ELISA Kit



Catalog No: #EK19858

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)

Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

## Description

Product Name	Human CHRNA2/Neuronal acetylcholine receptor subunit alpha-2 ELISA Kit
Brief Description	ELISA Kits
Applications	ELISA
Species Reactivity	Hu
Specificity	Natural and recombinant Human Neuronal acetylcholine receptor subunit alpha-2
Target Name	Neuronal acetylcholine receptor subunit alpha-2
Accession No.	Q15822
Uniprot	Q15822
GeneID	1135;

## Application Details

Assay Type: Sandwich Detect Range: 0.625-40 ng/mL Sensitivity: 0.270 ng/mL Sample Type: Serum, plasma, tissue homogenates and other biological fluids.

Sample Volume: 100ul Assay Time: 3-5h Detection wavelength: 450 nm gene\_full\_name: CHRNA2

Note: This product is for in vitro research use only