

Mouse Corin/Atrial natriuretic peptide-converting enzyme ELISA Kit



Catalog No: #EK19859

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Mouse Corin/Atrial natriuretic peptide-converting enzyme ELISA Kit
Brief Description	ELISA Kits
Applications	ELISA
Species Reactivity	Ms
Specificity	Natural and recombinant Mouse Atrial natriuretic peptide-converting enzyme
Target Name	Atrial natriuretic peptide-converting enzyme
Other Names	Corin,Low density lipoprotein receptor-related protein 4,Pro-ANP-converting enzyme,Crn,Lrp4
Accession No.	Q9Z319
Uniprot	Q9Z319
GeneID	53419;

Application Details

Assay Type: SandwichDetect Range:0.312-20ng/mL Sensitivity:0.183ng/mL Sample Type:Serum, plasma, tissue homogenates and other biological fluids.

Sample Volume:100ul Assay Time:3-5h Detection wavelength:450 nmgene_full_name:Corin

Note: This product is for in vitro research use only