

Human CTSS/Cathepsin S ELISA Kit

Catalog No: #EK19860



Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

Description

Product Name	Human CTSS/Cathepsin S ELISA Kit
Brief Description	ELISA Kits
Applications	ELISA
Species Reactivity	Hu
Specificity	Natural and recombinant Human Cathepsin S
Target Name	Cathepsin S
Accession No.	P25774
Uniprot	P25774
GeneID	1520;

Application Details

Assay Type: Sandwich Detect Range: 0.156-10ng/mL Sensitivity: 0.08ng/mL Sample Type: Serum, plasma, tissue homogenates and other biological fluids.
Sample Volume: 100ul Assay Time: 3-5h Detection wavelength: 450 nm gene_full_name: CTSS

Note: This product is for in vitro research use only