

Bovine AANAT/Serotonin N-acetyltransferase ELISA Kit



Catalog No: #EK19861

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Bovine AANAT/Serotonin N-acetyltransferase ELISA Kit
Brief Description	ELISA Kits
Applications	ELISA
Species Reactivity	Bovine
Specificity	Natural and recombinant Bovine Serotonin N-acetyltransferase
Target Name	Serotonin N-acetyltransferase
Other Names	Serotonin acetylase,Aralkylamine N-acetyltransferase,AA-NAT,SNAT
Accession No.	O02785
Uniprot	O02785
GeneID	281583;

Application Details

Assay Type: SandwichDetect Range:78-5000pg/mL Sensitivity:40.5pg/mL Sample Type:Serum, plasma, tissue homogenates and other biological fluids.
Sample Volume:100ul Assay Time:3-5h Detection wavelength:450 nmgene_full_name:AANAT

Note: This product is for in vitro research use only