

Mouse Aanat/Serotonin N-acetyltransferase ELISA Kit

Catalog No: #EK19862

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Mouse Aanat/Serotonin N-acetyltransferase ELISA Kit
Brief Description	ELISA Kits
Applications	ELISA
Species Reactivity	Ms
Specificity	Natural and recombinant Mouse Serotonin N-acetyltransferase
Target Name	Serotonin N-acetyltransferase
Other Names	Serotonin acetylase,Aralkylamine N-acetyltransferase,AA-NAT,Snat
Accession No.	O88816
Uniprot	O88816
GeneID	11298;

Application Details

Assay Type: Sandwich Detect Range:0.156-10ng/mL Sensitivity:0.082ng/mL Sample Type:Serum, plasma, tissue homogenates and other biological fluids.

Sample Volume:100ul Assay Time:3-5h Detection wavelength:450 nm gene_full_name:Aanat

Note: This product is for in vitro research use only