

Mouse Apex1/DNA-(apurinic or apyrimidinic site) lyase ELISA Kit



Catalog No: #EK19863

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Mouse Apex1/DNA-(apurinic or apyrimidinic site) lyase ELISA Kit
Brief Description	ELISA Kits
Applications	ELISA
Species Reactivity	Ms
Specificity	Natural and recombinant Mouse DNA-(apurinic or apyrimidinic site) lyase
Target Name	DNA-(apurinic or apyrimidinic site) lyase
Other Names	APEX nuclease, APEN, Apurinic-apyrimidinic endonuclease 1, AP endonuclease 1, REF-1, Redox factor-1, Ape, Apex, Ref1
Accession No.	P28352
Uniprot	P28352
GeneID	11792;

Application Details

Assay Type: Sandwich Detect Range: 0.156-10ng/mL Sensitivity: 0.083ng/mL Sample Type: Serum, plasma, tissue homogenates and other biological fluids.

Sample Volume: 100ul Assay Time: 3-5h Detection wavelength: 450 nm gene_full_name: Apex1

Note: This product is for in vitro research use only