



Mouse Cdkn2c/Cyclin-dependent kinase 4 inhibitor C ELISA Kit

Catalog No: #EK19864

Orders: order@signalwayantibody.com
Support: tech@signalwayantibody.com

Description

Product Name	Mouse Cdkn2c/Cyclin-dependent kinase 4 inhibitor C ELISA Kit
Brief Description	ELISA Kits
Applications	ELISA
Species Reactivity	Ms
Specificity	Natural and recombinant Mouse Cyclin-dependent kinase 4 inhibitor C
Target Name	Cyclin-dependent kinase 4 inhibitor C
Other Names	Cyclin-dependent kinase 6 inhibitor,p18-INK4c,p18-INK6
Accession No.	Q60772
Uniprot	Q60772
GeneID	12580;

Application Details

Assay Type: SandwichDetect Range:78-5000pg/mL Sensitivity:40.9pg/mL Sample Type:Serum, plasma, tissue homogenates and other biological fluids.
Sample Volume:100ul Assay Time:3-5h Detection wavelength:450 nmgene_full_name:Cdkn2c

Note: This product is for in vitro research use only