

Mouse Acr/Acrosin ELISA Kit

Catalog No: #EK19865

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Mouse Acr/Acrosin ELISA Kit
Brief Description	ELISA Kits
Applications	ELISA
Species Reactivity	Ms
Specificity	Natural and recombinant Mouse Acrosin
Target Name	Acrosin
Accession No.	P23578
Uniprot	P23578
GeneID	11434;

Application Details

Assay Type: Sandwich Detect Range: 78-5000pg/mL Sensitivity: 41pg/mL Sample Type: Serum, plasma, tissue homogenates and other biological fluids. Sample Volume: 100ul Assay Time: 3-5h Detection wavelength: 450 nm gene_full_name: Acr

Note: This product is for in vitro research use only