

## Rat Acr/Acrosin ELISA Kit

Catalog No: #EK19866

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

## Description

Product Name	Rat Acr/Acrosin ELISA Kit
Brief Description	ELISA Kits
Applications	ELISA
Species Reactivity	Rt
Specificity	Natural and recombinant Rat Acrosin
Target Name	Acrosin
Accession No.	P29293
Uniprot	P29293
GeneID	24163;

## Application Details

Assay Type: Sandwich Detect Range: 1.56-100ng/mL Sensitivity: 0.78ng/mL Sample Type: Serum, plasma, tissue homogenates and other biological fluids.

Sample Volume: 100ul Assay Time: 3-5h Detection wavelength: 450 nm gene\_full\_name: Acr

Note: This product is for in vitro research use only