

Mouse Ager/Advanced glycosylation end product-specific receptor ELISA Kit

Catalog No: #EK19867

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Mouse Ager/Advanced glycosylation end product-specific receptor ELISA Kit
Brief Description	ELISA Kits
Applications	ELISA
Species Reactivity	Ms
Specificity	Natural and recombinant Mouse Advanced glycosylation end product-specific receptor
Target Name	Advanced glycosylation end product-specific receptor
Other Names	Receptor for advanced glycosylation end products,Rage
Accession No.	Q62151
Uniprot	Q62151
GeneID	11596;

Application Details

Assay Type: Sandwich Detect Range:0.156-10ng/mL Sensitivity:0.088ng/mL Sample Type:Serum, plasma, tissue homogenates and other biological fluids.
Sample Volume:100ul Assay Time:3-5h Detection wavelength:450 nm gene_full_name:Ager

Note: This product is for in vitro research use only