

Rat Ager/Advanced glycosylation end product-specific receptor ELISA Kit



Catalog No: #EK19868

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Rat Ager/Advanced glycosylation end product-specific receptor ELISA Kit
Brief Description	ELISA Kits
Applications	ELISA
Species Reactivity	Rt
Specificity	Natural and recombinant Rat Advanced glycosylation end product-specific receptor
Target Name	Advanced glycosylation end product-specific receptor
Other Names	Receptor for advanced glycosylation end products,Rage
Accession No.	Q63495
Uniprot	Q63495

Application Details

Assay Type: Sandwich Detect Range:78-5000pg/mL Sensitivity:21pg/mL Sample Type:Serum, plasma, tissue homogenates and other biological fluids. Sample Volume:100ul Assay Time:3-5h Detection wavelength:450 nm gene_full_name:Ager

Note: This product is for in vitro research use only