

Mouse Cbs/Cystathionine beta-synthase ELISA Kit

Catalog No: #EK19869

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Mouse Cbs/Cystathionine beta-synthase ELISA Kit
Brief Description	ELISA Kits
Applications	ELISA
Species Reactivity	Ms
Specificity	Natural and recombinant Mouse Cystathionine beta-synthase
Target Name	Cystathionine beta-synthase
Other Names	Beta-thionase,Serine sulfhydrase
Accession No.	Q91WT9
Uniprot	Q91WT9
GeneID	12411;

Application Details

Assay Type: Sandwich Detect Range:0.156-10ng/mL Sensitivity:0.098ng/mL Sample Type:Serum, plasma, tissue homogenates and other biological fluids.

Sample Volume:100ul Assay Time:3-5h Detection wavelength:450 nm gene_full_name:Cbs

Note: This product is for in vitro research use only