

Mouse Cd34/Hematopoietic progenitor cell antigen CD34 ELISA Kit

Catalog No: #EK19870

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

| | |
|--------------------|--|
| Product Name | Mouse Cd34/Hematopoietic progenitor cell antigen CD34 ELISA Kit |
| Brief Description | ELISA Kits |
| Applications | ELISA |
| Species Reactivity | Ms |
| Specificity | Natural and recombinant Mouse Hematopoietic progenitor cell antigen CD34 |
| Target Name | Hematopoietic progenitor cell antigen CD34 |
| Other Names | CD34 |
| Accession No. | Q64314 |
| Uniprot | Q64314 |
| GeneID | 12490; |

Application Details

Assay Type: Sandwich Detect Range: 0.156-10ng/mL Sensitivity: 0.089ng/mL Sample Type: Serum, plasma, tissue homogenates and other biological fluids.

Sample Volume: 100ul Assay Time: 3-5h Detection wavelength: 450 nm gene_full_name: Cd34

Note: This product is for in vitro research use only