

Bovine ALOX12/Arachidonate 12-lipoxygenase, 12S-type ELISA Kit



Catalog No: #EK19871

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Bovine ALOX12/Arachidonate 12-lipoxygenase, 12S-type ELISA Kit
Brief Description	ELISA Kits
Applications	ELISA
Species Reactivity	Bovine
Specificity	Natural and recombinant Bovine Arachidonate 12-lipoxygenase, 12S-type
Target Name	Arachidonate 12-lipoxygenase, 12S-type
Other Names	12S-LOX, 12S-lipoxygenase
Accession No.	P27479
Uniprot	P27479
GeneID	282139;

Application Details

Assay Type: Sandwich Detect Range: 0.78-50ng/mL Sensitivity: 0.403ng/mL Sample Type: Serum, plasma, tissue homogenates and other biological fluids.
Sample Volume: 100ul Assay Time: 3-5h Detection wavelength: 450 nm gene_full_name: ALOX12

Note: This product is for in vitro research use only