

# Human ALOX12/Arachidonate 12-lipoxygenase, 12S-type ELISA Kit



Catalog No: #EK19872

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

## Description

Product Name	Human ALOX12/Arachidonate 12-lipoxygenase, 12S-type ELISA Kit
Brief Description	ELISA Kits
Applications	ELISA
Species Reactivity	Hu
Specificity	Natural and recombinant Human Arachidonate 12-lipoxygenase, 12S-type
Target Name	Arachidonate 12-lipoxygenase, 12S-type
Other Names	12S-LOX,Lipoxin synthase 12-LO,Platelet-type lipoxygenase 12,12LO,LOG12
Accession No.	P18054
Uniprot	P18054
GeneID	239;

## Application Details

Assay Type: Sandwich Detect Range:1.56-100ng/mL Sensitivity:0.78ng/mL Sample Type:Serum, plasma, tissue homogenates and other biological fluids.  
Sample Volume:100ul Assay Time:3-5h Detection wavelength:450 nm gene\_full\_name:ALOX12

Note: This product is for in vitro research use only