

Mouse Alox12/Arachidonate 12-lipoxygenase, 12S-type ELISA Kit



Catalog No: #EK19873

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Mouse Alox12/Arachidonate 12-lipoxygenase, 12S-type ELISA Kit
Brief Description	ELISA Kits
Applications	ELISA
Species Reactivity	Ms
Specificity	Natural and recombinant Mouse Arachidonate 12-lipoxygenase, 12S-type
Target Name	Arachidonate 12-lipoxygenase, 12S-type
Other Names	12S-LOX,Lipoxin synthase 12-LO,Platelet-type lipoxygenase 12,P-12LO,Alox12p
Accession No.	P39655
Uniprot	P39655
GeneID	11684;

Application Details

Assay Type: Sandwich Detect Range:0.156-10ng/ml Sensitivity:0.098ng/mL Sample Type:Serum, plasma, tissue homogenates and other biological fluids.
Sample Volume:100ul Assay Time:3-5h Detection wavelength:450 nm gene_full_name:Alox12

Note: This product is for in vitro research use only