

Mouse Apoa4/Apolipoprotein A-IV ELISA Kit

Catalog No: #EK19874

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Mouse Apoa4/Apolipoprotein A-IV ELISA Kit
Brief Description	ELISA Kits
Applications	ELISA
Species Reactivity	Ms
Specificity	Natural and recombinant Mouse Apolipoprotein A-IV
Target Name	Apolipoprotein A-IV
Other Names	Apo-AIV, Apolipoprotein A4
Accession No.	P06728
Uniprot	P06728
GeneID	11808;

Application Details

Assay Type: Sandwich Detect Range: 1.56-100ng/mL Sensitivity: 0.78ng/mL Sample Type: Serum, plasma, tissue homogenates and other biological fluids.
Sample Volume: 100ul Assay Time: 3-5h Detection wavelength: 450 nm gene_full_name: Apoa4

Note: This product is for in vitro research use only