

Rat ApoA4/Apolipoprotein A-IV ELISA Kit

Catalog No: #EK19875

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Rat ApoA4/Apolipoprotein A-IV ELISA Kit
Brief Description	ELISA Kits
Applications	ELISA
Species Reactivity	Rt
Specificity	Natural and recombinant Rat Apolipoprotein A-IV
Target Name	Apolipoprotein A-IV
Other Names	Apo-AIV, Apolipoprotein A4
Accession No.	P02651
Uniprot	P02651
GeneID	25080;

Application Details

Assay Type: Sandwich Detect Range: 0.312-20ng/mL Sensitivity: 0.171ng/mL Sample Type: Serum, plasma, tissue homogenates and other biological fluids.

Sample Volume: 100ul Assay Time: 3-5h Detection wavelength: 450 nm gene_full_name: ApoA4

Note: This product is for in vitro research use only