

Mouse C8a/Complement component C8 alpha chain ELISA Kit



Catalog No: #EK19877

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Mouse C8a/Complement component C8 alpha chain ELISA Kit
Brief Description	ELISA Kits
Applications	ELISA
Species Reactivity	Ms
Specificity	Natural and recombinant Mouse Complement component C8 alpha chain
Target Name	Complement component C8 alpha chain
Other Names	Complement component 8 subunit alpha
Accession No.	Q8K182
Uniprot	Q8K182
GeneID	230558;

Application Details

Assay Type: Sandwich Detect Range: 0.312-20ng/mL Sensitivity: 0.177ng/mL Sample Type: Serum, plasma, tissue homogenates and other biological fluids.

Sample Volume: 100ul Assay Time: 3-5h Detection wavelength: 450 nm gene_full_name: C8a

Note: This product is for in vitro research use only