

# Rat Adamts1/A disintegrin and metalloproteinase with thrombospondin motifs 1 ELISA Kit

Catalog No: #EK19878

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

## Description

|                    |  |
|--------------------|--|
| Product Name       | Rat Adamts1/A disintegrin and metalloproteinase with thrombospondin motifs 1 ELISA Kit       |
| Brief Description  | ELISA Kits   |
| Applications       | ELISA  |
| Species Reactivity | Rt   |
| Specificity        | Natural and recombinant Rat A disintegrin and metalloproteinase with thrombospondin motifs 1 |
| Target Name        | A disintegrin and metalloproteinase with thrombospondin motifs 1                             |
| Other Names        | ADAM-TS 1  |
| Accession No.      | Q9WUQ1   |
| Uniprot            | Q9WUQ1   |

## Application Details

Assay Type: Sandwich Detect Range: 0.312-20ng/mL Sensitivity: 0.18ng/mL Sample Type: Serum, plasma, tissue homogenates and other biological fluids.

Sample Volume: 100ul Assay Time: 3-5h Detection wavelength: 450 nm gene\_full\_name: Adamts1

Note: This product is for in vitro research use only