

# Mouse Adamts7/A disintegrin and metalloproteinase with thrombospondin motifs 7 ELISA Kit



Catalog No: #EK19879

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

## Description

Product Name	Mouse Adamts7/A disintegrin and metalloproteinase with thrombospondin motifs 7 ELISA Kit
Brief Description	ELISA Kits
Applications	ELISA
Species Reactivity	Ms
Specificity	Natural and recombinant Mouse A disintegrin and metalloproteinase with thrombospondin motifs 7
Target Name	A disintegrin and metalloproteinase with thrombospondin motifs 7
Other Names	ADAM-TS 7,COMPase
Accession No.	Q68SA9
Uniprot	Q68SA9
GeneID	108153;

## Application Details

Assay Type: Sandwich Detect Range:0.156-10ng/mL Sensitivity:0.078ng/mL Sample Type:Serum, plasma, tissue homogenates and other biological fluids.  
Sample Volume:100ul Assay Time:3-5h Detection wavelength:450 nm gene\_full\_name:Adamts7

Note: This product is for in vitro research use only