



Rat Adamts7/A disintegrin and metalloproteinase with thrombospondin motifs 7 ELISA Kit

Catalog No: #EK19880

Orders: order@signalwayantibody.com
Support: tech@signalwayantibody.com

Description

Product Name	Rat Adamts7/A disintegrin and metalloproteinase with thrombospondin motifs 7 ELISA Kit
Brief Description	ELISA Kits
Applications	ELISA
Species Reactivity	Rt
Specificity	Natural and recombinant Rat A disintegrin and metalloproteinase with thrombospondin motifs 7
Target Name	A disintegrin and metalloproteinase with thrombospondin motifs 7
Other Names	ADAM-TS 7,COMPase
Accession No.	Q1EHB3
Uniprot	Q1EHB3
GeneID	315879;

Application Details

Assay Type: SandwichDetect Range:0.312-20ng/mL Sensitivity:0.183ng/mL Sample Type:Serum, plasma, tissue homogenates and other biological fluids.
Sample Volume:100ul Assay Time:3-5h Detection wavelength:450 nmgene_full_name:Adamts7

Note: This product is for in vitro research use only