

Human CHKA/Choline kinase alpha ELISA Kit

Catalog No: #EK19881

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Human CHKA/Choline kinase alpha ELISA Kit
Brief Description	ELISA Kits
Applications	ELISA
Species Reactivity	Hu
Specificity	Natural and recombinant Human Choline kinase alpha
Target Name	Choline kinase alpha
Other Names	CK,CHETK-alpha,Ethanolamine kinase,EK,CHK,CKI
Accession No.	P35790
Uniprot	P35790
GeneID	1119;

Application Details

Assay Type: Sandwich Detect Range:0.312-20ng/mL Sensitivity:0.156ng/mL Sample Type:Serum, plasma, tissue homogenates and other biological fluids.

Sample Volume:100ul Assay Time:3-5h Detection wavelength:450 nm gene_full_name:CHKA

Note: This product is for in vitro research use only