

Mouse Chka/Choline kinase alpha ELISA Kit

Catalog No: #EK19882

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Mouse Chka/Choline kinase alpha ELISA Kit
Brief Description	ELISA Kits
Applications	ELISA
Species Reactivity	Ms
Specificity	Natural and recombinant Mouse Choline kinase alpha
Target Name	Choline kinase alpha
Other Names	CK, CHETK-alpha, Ethanolamine kinase, EK, Chk
Accession No.	O54804
Uniprot	O54804
GeneID	12660;

Application Details

Assay Type: Sandwich Detect Range: 0.156-10ng/mL Sensitivity: 0.08ng/mL Sample Type: Serum, plasma, tissue homogenates and other biological fluids.

Sample Volume: 100ul Assay Time: 3-5h Detection wavelength: 450 nm gene_full_name: Chka

Note: This product is for in vitro research use only