

Rat Chka/Choline kinase alpha ELISA Kit

Catalog No: #EK19883

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Rat Chka/Choline kinase alpha ELISA Kit
Brief Description	ELISA Kits
Applications	ELISA
Species Reactivity	Rt
Specificity	Natural and recombinant Rat Choline kinase alpha
Target Name	Choline kinase alpha
Other Names	CK,CHETK-alpha,Ethanolamine kinase,EK,Chk,Ckr
Accession No.	Q01134
Uniprot	Q01134
GeneID	29194;

Application Details

Assay Type: Sandwich Detect Range:0.78-50ng/ml Sensitivity:0.19ng/ml Sample Type:Serum, plasma, tissue homogenates and other biological fluids. Sample Volume:100ul Assay Time:3-5h Detection wavelength:450 nm gene_full_name:Chka

Note: This product is for in vitro research use only