

Bovine ARRB2/Beta-arrestin-2 ELISA Kit

Catalog No: #EK19886

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Bovine ARRB2/Beta-arrestin-2 ELISA Kit
Brief Description	ELISA Kits
Applications	ELISA
Species Reactivity	Bovine
Specificity	Natural and recombinant Bovine Beta-arrestin-2
Target Name	Beta-arrestin-2
Other Names	Arrestin beta-2, Arrestin-3
Accession No.	P32120
Uniprot	P32120

Application Details

Assay Type: Sandwich Detect Range: 0.312-20ng/mL Sensitivity: 0.156ng/mL Sample Type: Serum, plasma, tissue homogenates and other biological fluids.

Sample Volume: 100ul Assay Time: 3-5h Detection wavelength: 450 nm gene_full_name: ARRB2

Note: This product is for in vitro research use only